

Attendant Control Keypad - ACK **Data sheet**



ACK

With the OpenBus™ system it is possible to use ACK membrane front covers as attendant controls or hand controls integrated in the bed side rails.

There are two different variants of ACK: ACK1 and ACK3. The ACK1 is a single membrane front cover, whereas the ACK3 comes with two membrane front covers, typically used on an inside side rail and an outside side rail.

ACK offers easy access to critical positioning functions and correct bed position in case of emergency. Another benefit is high adaptation to customer needs with customised key and diode layouts, and membrane front covers adapted to the requested design and colour.

Further benefits are included in the ACK3 version, as this offers double and flexible bed control for the patient and the nurse to give increased comfort and efficiency. ACK3 allows for full intuitive bed control with up to 2x15 keys and 2x7 LEDs.



Features and options:


- The membrane front covers can be used as controls for nurse and patient
- Placed conveniently in the bed side rails allowing optimum ease of use
- Up to 15 keys and 7 LEDs possible with the ACK1
- Up to 2x15 keys and 2x7 LEDs possible with the ACK3
- Protection against single fault condition according to IEC60601-1 when connected to a LINAK OpenBus™ control box
- An LED which indicates when the battery is being charged (yellow)
- There is no IP rating as standard. The customer must ensure a proper application design and build-in of the ACK unit to achieve the correct IP degree (e.g. IPX4).
- Straight cables: 1250, 1800 or 2500 mm
- Connection to CB or MJB via a "modular jack" cable (6-wired)
- The standard ACK colour is grey (RAL 7035)
- Possible to customise design and colours of the membrane front covers
- Possible to customise the layout of the keys and LED's
- Possible to customise SW (LixEdit)

Usage:

- Operation temperature: +5 °C to +40 °C
- Storage temperature: -10 °C to +50 °C
- Compatibility: Compatible with LINAK control boxes. Please contact LINAK.
- Relative humidity: 20% to 80% non-condensing
- Atmospheric pressure: 700 to 1060 hPa
- Meters above sea level: Max. 3000 meters
- Approvals: IEC60601-1

ANSI/AAMI ES60601-1 CAN/CSA-22.2 No 60601-1

ACK Ordering example:



An electronic interface (PCBA), a membrane front cover and a cable are included when ordering an ACK – Attendant Control Keypad.

Overview


Please use the below membrane front cover and functionality guide for easy ordering

| ACK | Front cover no. | Functionality | Keys | LEDs | Graphics for ACK standard membrane front covers |
|------|------------------------------------|---------------|------|------|---|
| ACK1 | S0 = 10ACK00S28 | 04 | 11 | 5 | LINAK ES |
| ACK1 | S4 = 10ACK00S4 | 01 | 6 | 3 | LINAX 23 |
| ACK1 | S5 = 10ACK00S5 | 02 | 11 | 5 | |
| ACK1 | \$6 = 10ACK00\$6 | 03 | 6 | 3 | LRUNCES C |
| ACK3 | S7 = 10ACK00S5 + 10ACK00S4 | 21 | 11+6 | 5+3 | UINACES UINACES |
| ACK3 | S8 = 10ACK00S6 + 10ACK00S28 | 23 | 6+11 | 3+5 | LINAK 25 |
| ACK3 | SA = 10ACK00S96* + 10ACK00S4 | 21 | 11+6 | 5+3 | LIMAK 58 |
| ACK3 | SB = 10ACK00S6 + 10ACK00S97* | 23 | 6+11 | 3+5 | LINAK 28 |


^{*} with long tail (see dimensions)

Tail connectors

Key tail




LED tail



Dimensions


Patient membrane front cover



Keep a minimum Creepage of 2.5 mm to all internal wiring. (EN60601-1:2006)


Best practice! Wiring/Net nearest to Creepage border should be GND.


Drawing No. 0971003




Drawing No. 0971001

Nurse membrane front cover






| P/N table | | |
|------------|-----|-----|
| LINAK P/N | L1 | L2 |
| 10ACK00S28 | 92 | 72 |
| 10ACK00S97 | 157 | 137 |



Standard front covers


Front cover no.: 10ACK00S4




Front cover no.: 10ACK00S5 / 10ACK00S96



Front cover no.: 10ACK00S6



Front cover no.: 10ACK00S28 / 10ACK00S97



For customised membrane front covers please contact LINAK A/S.

Copyright © LINAK 2021.08 MA M9-02-216-H. Chapter 9.17

LINAK" accepts no responsibility for possible errors or inaccuracies in catalogues, brochures, and other material. LINAK reserves the right to change its products without prior notice. LINAK cannot guarantee product availability and reserves the right to discontinue the sale of any product. The user is responsible for determining the suitability of LINAK products for a specific application. All sales are subject to the 'Standard Terms of Sale and Delivery', available on LINAK websites.

LINAK and the LINAK logotype are registered trademarks of LINAK A/S. All rights reserved.