

Hand Control HB80

Data Sheet



Contents

Introduction.....	3
Features and options	4
Usage	4
Dimensions	5
HB81 - HB85 Dimensions (mm).....	5
HB86 dimensions (mm).....	5
Ordering example	6
Front covers	7
HB8x Standard front covers, transparent	8
HB8x Standard front covers, not transparent	8
HB8x Standard front covers, not transparent	9
HB8x Standard front covers, not transparent	10
HB8x Antimicrobial front cover, example.....	10
HB80 Standard front covers	11
Information	11
Contacts	12

Introduction



The HB80 hand control has an optimised ergonomic design shaped for the hand. The hand control is suitable for all kinds of MEDLINE® and CARELINE® applications, such as hospital beds, patient lifts, treatment and examination couches etc.

The HB80 hand control is available in versions with up to 10 or 12 activation buttons.

Further benefits are:

Optional antimicrobial HB80 material (housing, hook, front cover) to limit the growth of microbes.

Available in several different standard front cover versions with a variation of symbols. Customised front covers can easily be made for a differentiated look.

Optional blue or green background lighting for use with LINAK OpenBus™ control boxes. For differentiation and easy usage in dark surroundings.

High water resistance with optional IP degree of IPX6 Washable DURA™.

The HB80 is integrated in the LINAK Configurator Econ Smart Client.

Features and options

Hand control for up to 6 actuators

Housing colour:	Light grey or dark grey
Protection class:	IPX4, IPX6 Washable DURA™
Number of buttons:	Up to 10 or 12
Front cover:	Standard or customised
HB80 antimicrobial:	Yes, as an option
Cable:	Coiled or straight
Button type:	Easy touch tactile switches
Memory option:	Possible
Background lighting:	Blue or green (for use with OpenBus™ control)
Attach:	Via hook or magnet

Usage



Operation temperature:	+5 °C to +40 °C
Storage temperature:	-10 °C to +50 °C
Compatibility:	Compatible with many LINAK control boxes. Please contact LINAK.
Relative humidity:	20% to 80% – non-condensing
Atmospheric pressure:	700 to 1060 hPa
Meters above sea level:	Max. 3000 meters
Approvals:	IEC60601-1 ANSI/AAMI ES60601-1 CAN/CSA-22.2 No 60601-1

The HB86 version has a shorter distance between the buttons and cannot be approved according to EN IEC60601-2-52 Application Environment 4 for care beds used in Domestic area (or EN1970).

HB80 is designed and tested in accordance with EN60601-2-52 cl. 201.11.6.6.101 (machine washable medical beds).

Dimensions

HB81 - HB85 Dimensions (mm)





Hook option



Magnet option

HB86 dimensions (mm)



Hook option



Magnet option

Ordering example


[illegible]

The antimicrobial HB80 version includes active additives in the plastic of the hand control housing and the hook. The front cover has a second layer that is antimicrobial.


Front covers

HB8x Standard front covers, transparent


(only for use with background light)*




00HB821D-020




00HB831D-021




00HB831D-022




00HB841D-001
prev. 00HB841A-001




00HB841D-002
prev. 00HB841A-002



00HB851D-002




00HB851D-003




00HB851D-004




00HB851D-007
prev. 00HB841A002




00HB851D-011




00HB851D-019




00HB851D-021



00HB851D-022



00HB851D-035




00HB861B-001
prev. 00HB861B001


* If a transparent cover is wanted without a backlight, a special article must be created.

HB8x Standard front covers, transparent


(only for use with background light)




00HB861B-002
prev. 00HB861B002



00HB861B-005




00HB861B-007




00HB861B-009

HB8x Standard front covers, not transparent

(not for use with background light)*



00HB811A-001



00HB811A-004




00HB811A-033




00HB821A-001



00HB821A-002



00HB821A-008



00HB821A-041



00HB831A-001



00HB831A-002



00HB831A-009

* If a non-transparent cover is wanted with a backlight, a special article must be created.

HB8x Standard front covers, not transparent

(not for use with background light)*




00HB831A-011



00HB831A-029
For OpenBus
V0 for COL50



00HB841A-003



00HB841A-006




00HB841A-008




00HB841A-032



00HB851A-005




00HB851A-012



00HB851A-013



00HB851A-020




00HB851A-037




00HB851A-056



00HB851A-057



00HB851A-072




00HB861E-008

* If a non-transparent cover is wanted with a backlight, a special article must be created.

HB8x Standard front covers, not transparent

(not for use with background light)*



00HB861E-010
V2 hand control for
CBJC (phase-out is
ongoing for CBJC)

00HB861E-076

00HB861E-080

HB8x Antimicrobial front cover, example



00HB851N-002

HB80 Standard front covers



Information

When ordering check that the hand control functionality corresponds to the symbols on the front cover.

For all types attention should be given to ensure that the channels shown correspond to the channels on the chosen control box.

For customised front covers please contact LINAK A/S.

Contacts

FACTORIES

Denmark - Headquarters

LINAK A/S

Phone: +45 73 15 15 15

Fax: +45 74 45 80 48

Fax (Sales): +45 73 15 16 13

Web: www.linak.com

China

LINAK (Shenzhen) Actuator Systems, Ltd.

Phone: +86 755 8610 6656

Phone: +86 755 8610 6990

Web: www.linak.cn

Slovakia

LINAK Slovakia s.r.o.

Phone: +421 51 7563 444

Web: www.linak.sk

Thailand

LINAK APAC Ltd.

Phone: +66 33 265 400

Web: www.linak.com

USA

LINAK U.S. Inc.

Americas Headquarters

Phone: +1 502 253 5595

Fax: +1 502 253 5596

Web: www.linak-us.com
www.linak-latinamerica.com

SUBSIDIARIES

Australia

LINAK Australia Pty. Ltd

Phone: +61 3 8796 9777

Fax: +61 3 8796 9778

E-mail: sales@linak.com.au

Web: www.linak.com.au

Austria

LINAK Zweigniederlassung - Österreich (Wien)

Phone: +43 (1) 890 7446

Fax: +43 (1) 890 744615

E-mail: info@linak.de

Web: www.linak.at - www.linak.hu

Belgium

LINAK Actuator-Systems NV/SA

(Belgium & Luxembourg)

Phone: +32 (0)9 230 01 09

E-mail: beinfo@linak.be

Web: www.linak.be - www.fr.linak.be

Brazil

LINAK Do Brasil Comércio De Atuadores Ltda.

Phone: +55 (11) 2832 7070

Fax: +55 (11) 2832 7060

E-mail: info@linak.com.br

Web: www.linak.com.br

Canada

LINAK Canada Inc.

Phone: +1 502 253 5595

Fax: +1 416 255 7720

E-mail: info@linak.ca

Web: www.linak-us.com

Czech Republic

LINAK C&S s.r.o.

Phone: +42 058 174 1814

Fax: +42 058 170 2452

E-mail: info@linak.cz

Web: www.linak.cz - www.linak.sk

Denmark - International

LINAK International

Phone: +45 73 15 15 15

E-mail: info@linak.com

Web: www.linak.com

Denmark - Sales

LINAK Danmark A/S

Phone: +45 86 80 36 11

Fax: +45 86 82 90 51

E-mail: linak@linak-silkeborg.dk

Web: www.linak.dk

Finland

LINAK OY

Phone: +358 10 841 8700

E-mail: linak@linak.fi

Web: www.linak.fi

France

LINAK France E.U.R.L

Phone: +33 (0) 2 41 36 34 34

Fax: +33 (0) 2 41 36 35 00

E-mail: linak@linak.fr

Web: www.linak.fr

Germany

LINAK GmbH

Phone: +49 6043 9655 0

Fax: +49 6043 9655 60

E-mail: info@linak.de

Web: www.linak.de

India

LINAK A/S India Liaison Office

Phone: +91 120 4531797

Fax: +91 120 4786428

E-mail: info@linak.in

Web: www.linak.in

Italy

LINAK ITALIA S.r.l.

Phone: +39 02 48 46 33 66

Fax: +39 02 48 46 82 52

E-mail: info@linak.it

Web: www.linak.it

Japan

LINAK K.K.

Phone: 81-45-533-0802

Fax: 81-45-533-0803

E-mail: linak@linak.jp

Web: www.linak.jp

Malaysia

LINAK Actuators Sdn. Bhd.

Phone: +60 4 210 6500

Fax: +60 4 226 8901

E-mail: info@linak-asia.com

Web: www.linak.my

Netherlands

LINAK Actuator-Systems B.V.

Phone: +31 76 5 42 44 40 /

+31 76 200 11 10

E-mail: info@linak.nl

Web: www.linak.nl

New Zealand

LINAK New Zealand Ltd

Phone: +64 9580 2071

Fax: +64 9580 2072

E-mail: nzsales@linak.com.au

Web: www.linak.com.au

Norway

LINAK Norge AS

Phone: +47 32 82 90 90

E-mail: info@linak.no

Web: www.linak.no

Poland

LINAK Polska

LINAK Danmark A/S (Spółka Akcyjna)

Phone: +48 22 295 09 70 /

+48 22 295 09 71

E-mail: info@linak.pl

Web: www.linak.pl

Republic of Korea

LINAK Korea Ltd.

Phone: +82 2 6231 1515

Fax: +82 2 6231 1516

E-mail: info@linak.kr

Web: www.linak.kr

Slovakia

LINAK Slovakia S.R.O.

Phone: +421 51 7563 444

Web: www.linak.sk

Spain

LINAK Actuadores, S.L.u

Phone: +34 93 588 27 77

Fax: +34 93 588 27 85

E-mail: esma@linak.es

Web: www.linak.es

Sweden

LINAK Scandinavia AB

Phone: +46 8 732 20 00

Fax: +46 8 732 20 50

E-mail: info@linak.se

Web: www.linak.se

Switzerland

LINAK AG

Phone: +41 43 388 31 88

Fax: +41 43 388 31 87

E-mail: info@linak.ch

Web: www.linak.ch - www.fr.linak.ch

www.it.linak.ch

Taiwan

LINAK (Shenzhen) Actuator systems Ltd.

Taiwan Representative office

Phone: +886 2 272 90068

Fax: +886 2 272 90096

E-mail: sales@linak.com.tw

Web: www.linak.com.tw

Turkey

LINAK İth. İhr. San. ve Tic. A.Ş.

Phone: +90 312 4726338

Fax: +90 312 4726635

E-mail: info@linak.com.tr

Web: www.linak.com.tr

United Kingdom & Ireland

LINAK UK Limited

Phone: +44 (0)121 544 2211

Fax: +44 (0)121 544 2552

E-mail: sales@linak.co.uk

Web: www.linak.co.uk

DISTRIBUTORS

Argentina

Novotec Argentina SRL

Phone: 011-4303-8989 / 8900

Fax: 011-4032-0184

E-mail: info@novotecargentina.com

Web: www.novotecargentina.com

Colombia

MEM Ltda

Phone: +[57] (1) 334-7666

Fax: +[57] (1) 282-1684

E-mail: servicioalcliente@memltda.com.co

Web: www.mem.net.co

India

Mechatronics Control Equipments India Pvt Ltd

Phone: +91-44-28558484, 85

E-mail: bala@mechatronicscontrol.com

Web: www.mechatronicscontrol.com

Indonesia

PT. Himalaya Everest Jaya

Phone: +6 221 544 8956 /+6 221 544 8965

Fax: +6 221 619 1925

Fax (Sales): +6 221 619 4658

E-mail: hejplastic-div@centrin.net.id

Web: www.hej.co.id

Israel

NetivTech LTD

Phone: +972 55-2266-535

Fax: +972 2-9900-560

Email: info@NetivTech.com

Web: www.netivtech.com

Singapore

Servo Dynamics Pte Ltd

Phone: +65 6844 0288

Fax: +65 6844 0070

E-mail: servodynamics@servo.com.sg

South Africa

Industrial Specialised Applications CC

Phone: +27 011 466 0346

E-mail: gartht@isagroup.co.za

Web: www.isaza.co.za

United Arab Emirates

Mechatronics

Phone: +971 4 267 4311

Fax: +971 4 267 4312

E-mail: mechtron@emirates.net.ae

Terms of use

LINAK® takes great care in providing accurate and up-to-date information on its products. However, the user is responsible for determining the suitability of LINAK products for a specific application. Due to continual development, LINAK products are subject to frequent modifications and changes. LINAK reserves the rights to conduct modifications, updates, and changes without any prior notice. For the same reason, LINAK cannot guarantee the correctness and actual status of imprinted information on its products.

LINAK uses its best efforts to fulfil orders. However, for the reasons mentioned above, LINAK cannot guarantee availability of any particular product at any given time. LINAK reserves the right to discontinue the sale of any product displayed on its website or listed in its catalogues or in other written material created and produced by LINAK, LINAK subsidiaries, or LINAK affiliates. All sales are subject to the 'Standard Terms of Sale and Delivery for LINAK A/S' available on LINAK websites. LINAK and the LINAK logotype are registered trademarks of LINAK A/S. All rights reserved.